Detection Test for Staphylococcal Penicillinase Activity**Store at -20°C****Expiry date:** see packaging**40 tests****IVD****CE****EN**molsid.com/techdoc**REF RST1****1. INTENDED PURPOSE AND USE**

The ResiScreen® Staph product is a visual, qualitative *in vitro* diagnostic test intended for detection of β -lactamase enzymatic activity from *Staphylococcus* spp. bacteria.

It is intended to assist in the presumptive diagnosis of penicillin-resistant *Staphylococcus* species by providing physicians with additional information on β -lactamase enzyme activity in the sample. It serves as a complement to routine antimicrobial susceptibility testing and is not intended to replace the institution's established protocols.

The ResiScreen® Staph test is performed on pure, isolated colonies previously identified as a *Staphylococcus* spp. strain obtained after primary isolation from culture of patient sample.

The product is intended to be used by a laboratory professional, in healthcare settings. It is not automated, nor a companion diagnostic.

2. SUMMARY AND EXPLANATION

The ability of certain bacteria to produce enzymes that inactivate β -lactam antibiotics, i.e., penicillins and cephalosporins, by breaking down β -lactam ring, has been recognized since the 1940s with first discovery by Abraham and Chain (1). Today a variety of these β -lactamase enzymes with different substrate specificities have been isolated across numerous bacterial species and they constitute a major and clinically relevant mechanism of resistance to β -lactam antimicrobials (2,3). In *Staphylococcus aureus* and coagulase-negative Staphylococci species, resistant strains produce a single class of β -lactamase known as BlaZ. This enzyme selectively hydrolyses penicillinase-labile penicillins and represent a major contributor of β -lactam resistance in these bacteria (4,5).

Detecting β -lactamase activity in *Staphylococcus* species is key for identifying antibiotic-resistant strains, especially for methicillin-sensitive bacteria. This information can assist in guiding appropriate treatment decisions, ensuring the right antibiotics are used, and helping to prevent the spread of antimicrobial resistance. Various methods exist to detect β -lactamase activity, including chromogenic, iodometric and acidimetric assays.

3. PRINCIPLE

The ResiScreen® Staph test utilizes fluorescence as a readout of β -lactamase enzymatic activity in a *Staphylococcus* spp. sample. It is based on the chemical SmartID® technology biotracer made of 3 components: a fluorophore (ELF-97), a spacer and a β -lactamase substrate unit.

Upon contact with the targeted β -lactamase enzyme, the SmartID® probe is cleaved at a specific site on the substrate part leading to probe self-collapse and release of the fluorophore. Free ELF-97 precipitates, acquiring green fluorescence properties under UV light. Detection of fluorescence in a bacterial sample containing the probe indicates presence of the specific BlaZ β -lactamase enzymatic activity. In absence of active β -lactamase, the biotracer is not degraded and little to no fluorescence is emitted within the time limits of the test.

Detection Test for Staphylococcal Penicillinase Activity

4. REAGENTS

Ready to use Molsid SmartID® biotracer solution. One tube contains enough solution for 10 tests.

5. STORAGE

Store tubes in their original container at -20°C until use. After removing the required amount of solution for the tests, return the tube to -20°C. Keep tubes away from sources of heat.

In these conditions, ResiScreen® Staph solution is stable until at least the expiry date shown on the label.

In-use stability has been demonstrated for up to five thawing cycles of a single tube. For additional cycles, appropriate quality control should be performed in parallel; refer to the quality control section.

6. PRECAUTIONS AND WARNINGS

For *in vitro* diagnostic use. For professional use by laboratory personnel.

Observe aseptic methods and established precautions against microbiological hazards throughout all procedures.

After use, dispose of used samples, petri dishes and other contaminated materials according to local regulation and safety methods. The kit must be treated as biohazardous and chemically contaminated waste. Dispose of all materials in accordance with local, regional, and institutional regulations governing laboratory waste.

Refer to the Material Safety Data Sheet (MSDS) for information on potentially hazardous components.

Users should report any serious incident related to the device to the Manufacturer and National Competent Authority.

- Use Appropriate UV protective equipment when performing test revelation under UV light.
- Directions should be read and followed carefully. Deviations from the recommended protocol may result in false positive (*i.e.* from improper storage) or false negative results (*i.e.* from insufficient colony suspension in the mix with the biotracer, or inadequate UV lamp/system).
- This product should not be used if there is visible damage to the tube containing the biotracer or any other signs of deterioration.
- Do not use the product if the solution in the tube presents solid particles once it has been properly thawed.
- Do not use the product after expiry date.
- Product ingestion, inhalation, or contact with the skin or eyes should be avoided.

7. MATERIALS REQUIRED BUT NOT PROVIDED


(1) Pipette, (2) sterile pipette tips, (3) sterile microtubes, (4) sterile NaCl 0.85-0.90% (saline) solution, (5) Mueller Hinton agar plate, (6) sterile 1µL size inoculating loops, (7) UV_{365nm}-light transmitting lamp/table, (8) vortex, (9) incubator, (10) suitable (personal) protective equipment.

UV equipment shall be a **365nm light source** with a **power between 5 W and 50 W portable lamp or UV table**. The use of a UV darkroom is not recommended, as it may lead to false positive results during interpretation.

Detection Test for Staphylococcal Penicillinase Activity

8. PROCEDURE

Overview of the procedure

Test Procedure

The procedure is given for one (1) assay (control drop and test drop)

1. Remove the ResiScreen® Staph tube(s) from storage and allow the biotracer solution to thaw entirely at room temperature.
2. Mix by vortexing the tube.
3. Add **10µL** of biotracer solution into an empty microtube.
4. Prepare the bacterial suspension:
 - a. Take one **full** calibrated loop (**loop size 1µL**) of bacterial colonies (see Figure1) and resuspend it in **100µL** of 0.85%-0.90% saline solution. Use multiple colonies to avoid atypical variant.
 - b. Vortex and mix.


Figure1: Correct filling of the 1µl inoculating loop. Panel A: Partial loop filling with insufficient inoculum, may result in weak or absent fluorescence. Panel B: Complete loop filling with appropriate inoculum, required for valid test performance.

5. Add **10µL** of bacterial suspension in the tube containing the biotracer.
6. Vortex for 5 seconds.
7. Gently drop **10µL** of bacterial suspension directly on a **Mueller Hinton agar plate** (control drop).
8. Gently drop **10µL** of the mix bacterial suspension + biotracer on the same agar plate (test drop). Place both drop next to each other but ensuring no contact to avoid any contamination.
9. Wait for the drops to completely dry for roughly 5-10 minutes (do not reverse the plate immediately to prevent drop spillage or mixing).
10. Incubate the plate reversed for 16-24h (overnight) at 35°C ± 2 °C.
11. Read the plates:
 - a. Put the plate under a 365nm UV light source preferably on a dark surface (avoid white surface) or on a 365nm UV table.
 - b. Read both the control and the test drops. Interpret the result using the table in the interpretation section.

Detection Test for Staphylococcal Penicillinase Activity

9. INTERPRETATION

It is recommended to initially test strains specified in the Quality Control section to obtain an initial point of reference under user specific operating conditions.

The test is interpreted by looking for a bright green fluorescence on the test drop. The control drop is used to show potential growth issues.

Under the interpretation table are given some examples of what to visually expect from a test drop for a positive and a negative strain.

Control drop	Test drop	Result	Interpretation
No fluorescence	No fluorescence	Negative	Absence of β -lactamase activity
No fluorescence	Fluorescence	Positive	Presence of a β -lactamase activity
Fluorescence	No fluorescence	Uninterpretable*	
Fluorescence	Fluorescence		

*If the control drop shows fluorescence, the test is invalid, and retest must be considered.

Visual example of positive and negative results under UV_{365nm}:

Example of positive results:

Example of negative results:

(For each picture, the control drop is on the left and the test drop is on the right)

Variation between different references of Muller Hinton medium and Petri dish may influence the perceived fluorescence intensity.

For some species of coagulase-negative *Staphylococcus*, slower growth was observed in the presence of the biotracer with colonies instead of a bacterial carpet after 16-24h of incubation. It does not affect test results. Interpret fluorescence as you would do for a carpet.

For few coagulase-negative *Staphylococcus* strains, complete growth inhibition was observed. In this case, the test should be considered as uninterpretable and not considered for presumptive diagnosis of β -lactamase activity.

10. PERFORMANCE CHARACTERISTICS.

In an independent retrospective clinical study with 301 strains representative of the genetic diversity of staphylococcal clinical isolates (200 *S. aureus*, 101 coagulase-negative Staphylococci) there was 100% concordance between ResiScreen® Staph results and *blaZ* gene detection by whole genome sequencing method (*Clinical Performance report, VAR-P2302*). Of the 200 *S. aureus* isolate, 122 were *blaZ* positive and 78 were *blaZ* negative. Of the 101 coagulase negative staphylococci, 54 were *blaZ* positive and 47 were *blaZ* negative.

From this study, ResiScreen Staph® calculated performance results are as follows:

S. aureus (n=200)		
Parameter	Performance	CI95 range
Sensitivity	100%	[97%-100%]
Specificity	100%	[95.4%-100%]

Detection Test for Staphylococcal Penicillinase Activity

Concordance	100%	[98.2%-100%]
Coagulase negative staphylococci (n=101)		
Parameter	Performance	CI95 range
Sensitivity	100%	[93.4%-100%]
Specificity	100%	[92.5%-100%]
Concordance	100%	[96.4%-100%]

11. LIMITATIONS

ResiScreen® Staph test is designed to determine β-lactamase activity in *Staphylococcus* species as a complement to conventional antimicrobial susceptibility testing and should not be used as a replacement to antimicrobial susceptibility testing.

ResiScreen® Staph test is not to be used directly with human samples, with other sources containing mixed microbial flora or with bacteria other than *Staphylococcus* spp.

12. USER QUALITY CONTROL

Control reference cultures can be tested with each unknown sample using the following recommended strains. Systematic quality control of product functionality is not required for each run.

Test Strain	Expected Results
<i>Staphylococcus aureus</i> ATCC® 25923	Negative
<i>Staphylococcus aureus</i> ATCC® 29213	Positive

Quality control requirements must be performed in accordance with applicable regulations or accreditation requirements and your laboratory’s standard Quality Control procedures.

13. PACKAGING

4 tubes of ResiScreen® Staph.....10 tests/tube

14. REFERENCES

1. Abraham EP, Chain E. An enzyme from bacteria able to destroy penicillin. 1940. Rev Infect Dis. 1988 Jul-Aug;10(4):677-8. PMID: 3055168.
2. Bush K. Past and Present Perspectives on β-Lactamases. Antimicrob Agents Chemother. oct 2018;62(10):e01076-18.
3. Bush K, Bradford PA. Epidemiology of β-Lactamase-Producing Pathogens. Clin Microbiol Rev. 18 mars 2020;33(2):e00047-19.
4. Adams AP, Barry AL, Benner EJ. A Simple, Rapid Test to Differentiate Penicillin-susceptible from Penicillin-resistant *Staphylococcus aureus*.
5. Kirby WMM. Extraction of a Highly Potent Penicillin Inactivator from Penicillin Resistant *Staphylococci*. Science. 2 juin 1944;99(2579):452-3.

Detection Test for Staphylococcal Penicillinase Activity


Change History

Date	Version	Change Summary
2026-03	1.0	Initial release

Detection Test for Staphylococcal Penicillinase Activity

SYMBOL GLOSSARY

Some symbols listed below may not apply to this product.

-  Manufacturer
-  Authorized representative in Switzerland
-  Catalog number
-  Batch code
-  Use-by date
-  Storage Temperature limit
-  Consult electronic instructions for use.
-  Do not use if package is damaged.
-  *In vitro* diagnostic medical device
-  Contains sufficient amount for <n> tests.
-  Unique Device Identifier

Molsid SAS
60F avenue Rockefeller
69008 Lyon - FRANCE
Website: <https://molsid.com>
E-mail: contact@molsid.com

Casus Switzerland GmbH
Hinterbergstrasse 49
6312 Steinhausen
SWITZERLAND
E-mail: swissrep@casusconsulting.com

ATCC® is a trademark of American Type Culture Collection.
Molsid, Molsid logo, SmartID® and ResiScreen® are trademarks of Molsid.
All other trademarks are the property of their respective owners. © 2026 Molsid. All rights reserved.

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