Detection Test for Staphylococcal Penicillinase Activity

Store at -20°C

Expiry date: see packaging

ENGLISH

40 tests

IVD
CE

molsid.com/techdoc
REF
RST1

1. INTENDED USE

The ResiScreen® Staph product is a visual, qualitative test for detection of β -lactamase enzymatic activity from *Staphylococcus* spp. bacteria.

The ResiScreen® Staph test is performed on pure, isolated colonies previously identified as a *Staphylococcus* spp. strain obtained after primary isolation from culture of sample.

The test is for use by trained laboratory professionals, in healthcare settings. It is not automated, nor a companion diagnostic.

The result of ResiScreen® Staph is intended to give complementary information in addition to established antimicrobial susceptibility testing (AST) and to other relevant microbiological information.

ResiScreen® Staph is not intended to replace conventional AST.

2. SUMMARY AND EXPLANATION

The ability of certain bacteria to produce enzymes that inactivate β -lactam antibiotics, i.e., penicillins and cephalosporins, by breaking down β -lactam ring, has been recognized since the 1940s with first discovery by Abraham and Chain (1). Today a variety of these β -lactamase enzymes with different substrate specificities have been isolated across numerous bacterial species and they constitute a major and clinically relevant mechanism of resistance to β -lactam antimicrobials (2,3). In *Staphylococcus aureus* and coagulase-negative Staphylococci species, some strains produce a single class of β -lactamase known as BlaZ. This enzyme selectively hydrolyses penicillinase-labile penicillins and represent one of the molecular mechanisms that can contribute to penicillins resistance in these bacteria (4,5). However, test based on the detection of this enzyme cannot replace conventional susceptibility testing for β -lactams as other resistance mechanisms also influence the results of such tests, including *mecA/C* mediated methicillin resistance. Various methods, exist to detect β -lactamase activity, including iodometric, acidimetric or chromogenic assays such as nitrocefin based tests.

3. PRINCIPLE OF THE METHOD

The ResiScreen® Staph test utilizes fluorescence as a readout of a β -lactamase enzymatic activity produced by a *Staphylococcus* spp. strain. It is based on Molsid's SmartID® fluorogenic reagent 352600-C (the SmartID® biotracer) composed of the ELF-97 fluorophore linked to a cephalosporin-based substrate unit, recognised by β -lactamases.

Upon contact with a β -lactamase enzyme, the SmartID® biotracer is cleaved at a specific site leading to the release of the fluorophore. Free ELF-97 precipitates, acquiring green fluorescence properties under UV_{365nm} light. Detection of fluorescence in a bacterial sample containing the biotracer indicates the specific presence of a β -lactamase enzymatic activity. In absence of active β -lactamase, the biotracer is not degraded and no fluorescence is emitted within the time limits of the test.

4. REAGENTS

Ready to use Molsid SmartID® biotracer solution. One tube contains enough solution (115µl) for 10 tests.

Detection Test for Staphylococcal Penicillinase Activity

5. STORAGE

Store tubes in their original container at -20°C until use. **After removing the required amount of solution for the tests, return the tube to -20°C.** Keep tubes away from sources of heat.

In these conditions, ResiScreen® Staph solution is stable until at least the expiry date shown on the label.

In-use stability has been demonstrated for up to five thawing cycles of a single tube. For additional cycles, appropriate quality control should be performed in parallel of the unknown sample; refer to the quality control section.

6. PRECAUTIONS AND WARNINGS

For *in vitro* diagnostic use. For professional use by laboratory personnel.

Observe aseptic methods and established precautions against microbiological hazards throughout all procedures.

Dispose of all used samples, used reagents, and any other contaminated disposable materials according to local regulation and safety methods. It is the responsibility of each laboratory to handle solid and liquid waste according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

Refer to the Material Safety Data Sheet (MSDS) for information and for safe handling and disposal of the product.

Use appropriate UV protective equipment when performing test revelation under UV light.

Directions should be read and followed carefully. Deviations from the recommended protocol may result in false positive (e.g. from improper storage) or false negative results (e.g. from insufficient colony suspension in the mix with the biotracer, or inadequate UV lamp/system).

This product should not be used if there is visible damage to the tube containing the biotracer or any other signs of deterioration.

Do not use the product if the solution in the tube presents solid particles once it has been properly thawed.

Do not use the product after expiry date.

Product ingestion, inhalation, or contact with the skin or eyes should be avoided.

7. MATERIALS REQUIRED BUT NOT PROVIDED


(1) Pipette, (2) sterile pipette tips, (3) sterile microtubes, (4) sterile NaCl 0.85-0.90% (saline) solution, (5) Mueller Hinton agar plate, (6) sterile 1 µL size inoculating loops, (7) **UV_{365nm}-light transmitting lamp/table**, (8) vortex, (9) incubator, (10) suitable (personal) protective equipment.

⚠ UV equipment shall be a **365nm light source, such as a portable lamp or UV table**, with a **power between 5 W and 50 W**. The use of a UV darkroom is not recommended, as it may lead to false positive results during interpretation.

Detection Test for Staphylococcal Penicillinase Activity

8. PROCEDURE


Overview of the procedure

Test Procedure

The procedure is given for one (1) assay (control drop and test drop)

1. Remove the ResiScreen® Staph tube from storage and allow the biotracer solution to thaw entirely at room temperature.
2. Mix by vortexing the tube.
3. Add **10µL** of biotracer solution into an empty microtube.
4. Prepare the bacterial suspension:
 - a. Take one **full** calibrated loop (**loop size 1µL**) of bacterial colonies (see Figure1) and resuspend it in **100µL** of 0.85%-0.90% saline solution. Use multiple colonies to avoid atypical variant.
 - b. Vortex and mix.

*Figure1: Correct filling of the 1µl inoculating loop.
Panel A: Partial loop filling with insufficient inoculum, may result in weak or absent fluorescence.
Panel B: Complete loop filling with appropriate inoculum, required for valid test performance.*

5. Add **10µL** of bacterial suspension in the tube containing the biotracer.
6. Vortex for 5 seconds.
7. Gently drop **10µL** of the bacterial suspension - biotracer mixture onto a Mueller-Hinton agar plate. This is the **test drop**.
8. For the control drop: Gently drop **10µL** of the initial bacterial suspension (without biotracer) directly onto the same agar plate. Place both drop next to each other, while ensuring there is no contact to avoid any contamination.
9. Wait for the drops to completely dry for roughly 5-10 minutes (do not reverse the plate immediately to prevent drop spillage or mixing).
10. Incubate the plate reversed for 16-24h (overnight) at $35^{\circ}\text{C} \pm 2^{\circ}\text{C}$.

Detection Test for Staphylococcal Penicillinase Activity

11. Read the plates:

- a. Illuminate the bacterial drops directly with a 365 nm UV light source, ensuring the light reaches the drops directly, without passing through the culture medium.
 - **Torch-type UV light source:** Remove the lid and place the plate on a dark, non-reflective surface (do not use white surfaces).
 - **UV table lamp:** Remove the lid and position the plate upside down on the UV table surface, with the bacterial drops facing the UV light source.
- b. Read both the control and the test drops. Interpret the result using the table in the interpretation section.

9. INTERPRETATION

It is recommended to initially test strains specified in the Quality Control section to obtain an initial point of reference under user specific operating conditions.

The test is read by looking for a **bright green fluorescence** on the test drop under UV_{365nm}. The control drop is used to show potential growth issues.


Interpretation for the β -lactamase activity is given in the table below.

Table 1: Test interpretation

Control drop	Test drop	Result	Interpretation
No fluorescence	No fluorescence	Negative	Absence of β -lactamase activity
No fluorescence	Fluorescence	Positive	Presence of a β -lactamase activity
Fluorescence	No fluorescence	Uninterpretable*	
Fluorescence	Fluorescence		

*If the control drop shows fluorescence, the test is invalid, and a retest must be considered.

Figure 2: Example of visually expect positive and negative results under UV_{365nm}

Variation between different Muller Hinton medium and Petri dish references may influence the fluorescence intensity but does not affect the test result: any fluorescence, whether from a confluent lawn or from the border of the inoculation drop area, should be interpreted as a positive result.

For some *Staphylococcus* isolates, particularly coagulase-negative staphylococci (CoNS) (less than 3% of the tested panel), delayed growth may be observed after 16–24 hours of incubation in the presence of the biotracer. This results in a ring-shaped growth pattern or discrete colonies, rather than a confluent bacterial lawn covering the entire inoculation area. These growth patterns do not affect the validity of the result. Fluorescence should be interpreted as it would be for a confluent lawn (see fig. 2).

For a small number of CoNS strains, complete inhibition of bacterial growth may be observed in the presence of the biotracer (less than 1.3% of the tested panel). In this case, the result is invalid and cannot be used for β -lactamase activity detection.

Detection Test for Staphylococcal Penicillinase Activity

10. PERFORMANCE CHARACTERISTICS.

ResiScreen® Staph was evaluated independently in a retrospective clinical study. Three hundred and one (301) strains representative of the genetic diversity of staphylococcal clinical isolates (200 *S. aureus*, 101 coagulase-negative Staphylococci) were tested and compared to *blaZ* gene detection by whole genome sequencing - WGS, method (*Clinical Performance report, VAR_P2302_01*). Of the 200 *S. aureus* isolate, 122 were *blaZ* positive and 78 were *blaZ* negative by WGS. Of the 101 coagulase negative staphylococci, 54 were *blaZ* positive and 47 were *blaZ* negative by WGS.

In agreement with their *blaZ* status as determined by WGS, all isolates were correctly identified by the ResiScreen® Staph test. Study results are presented in table 1 below.

Table 2: Retrospective clinical performance study result.

Parameter	<i>S. aureus</i> (n=200)		Coagulase negative staphylococci (n=101)	
	Performance	CI95 range	Performance	CI95 range
Sensitivity	100%	[97%-100%]	100%	[93.4%-100%]
Specificity	100%	[95.4%-100%]	100%	[92.5%-100%]
Accuracy	100%	[98.2%-100%]	100%	[96.4%-100%]

11. LIMITATIONS

ResiScreen® Staph test is designed to provide an indication of β -lactamase enzymatic activity in a bacterial strain. Its performances have been validated on *Staphylococcus* species only and has not been evaluated on other bacterial species.

ResiScreen® Staph is to be performed on pure, isolated colonies previously identified as a *Staphylococcus* spp.. It is not intended for direct use on primary clinical specimens or with other sources or samples containing mixed microbial flora.

Detection of β -lactamase activity provides only supplementary information to AST on a specific β -lactam degradation mechanism (penicillinase) in *Staphylococcus* spp. and does not cover all resistance mechanisms to beta-lactams. The test result should not be interpreted as predictive of susceptibility to β -lactam antimicrobials and shall not replace comprehensive conventional antimicrobial susceptibility testing. ResiScreen® Staph is not intended for the detection or identification of *Staphylococcus* spp., including methicillin-resistant *Staphylococcus aureus* (MRSA).

ResiScreen® Staph result must not be used alone for establishing, modifying or discontinuing antimicrobial treatment. Conventional antimicrobial susceptibility testing and all relevant clinical and microbiological information must be obtained and considered as the primary basis for patient management.

12. USER QUALITY CONTROL

Systematic quality control of product functionality is not required for each run. When required, the following positive and negative reference strains are recommended as control strains.

Test Strain	Expected Results
<i>Staphylococcus aureus</i> ATCC® 25923	Negative
<i>Staphylococcus aureus</i> ATCC® 29213	Positive

Quality control requirements must be carried out in accordance with applicable regulations, accreditation standards and laboratory's Quality Control procedures.

Detection Test for Staphylococcal Penicillinase Activity

13. PACKAGING

4 tubes of ResiScreen® Staph.....10 tests/tube

14. REFERENCES

1. Abraham EP, Chain E. An enzyme from bacteria able to destroy penicillin. 1940. Rev Infect Dis. 1988 Jul-Aug;10(4):677-8. PMID: 3055168.
2. Bush K. Past and Present Perspectives on β -Lactamases. Antimicrob Agents Chemother. oct 2018;62(10):e01076-18.
3. Bush K, Bradford PA. Epidemiology of β -Lactamase-Producing Pathogens. Clin Microbiol Rev. 18 mars 2020;33(2):e00047-19.
4. Adams AP, Barry AL, Benner EJ. A Simple, Rapid Test to Differentiate Penicillin-susceptible from Penicillin-resistant Staphylococcus aureus.
5. Kirby WMM. Extraction of a Highly Potent Penicillin Inactivator from Penicillin Resistant Staphylococci. Science. 2 juin 1944;99(2579):452-3.

Technical Service and Support: contact your local Molsid's product distributor.

EU Only: Users shall report any serious incident related to the device to the Manufacturer (QRA@molsid.com) and National Competent Authority.

Change History		
Date	Version	Change Summary
2026-03	1.0	Initial release
2026-09	1.1	Editorial revision to improve clarity and precision of the sections: "intended use", "limitations", "test procedure" and "interpretation". Add contact for serious incident reporting for EU.

Detection Test for Staphylococcal Penicillinase Activity

SYMBOL GLOSSARY

Some symbols listed below may not apply to this product.

	Manufacturer
	Catalog number
	Batch code
	Use-by date
	Storage Temperature limit
	Consult electronic instructions for use.
	Do not use if package is damaged.
	<i>In vitro</i> diagnostic medical device
	Contains sufficient amount for <n> tests.
	Unique Device Identifier

Molsid SAS

60F avenue Rockefeller
69008 Lyon - FRANCE

Website: <https://molsid.com>

E-mail: contact@molsid.com

ATCC® is a trademark of American Type Culture Collection.

Molsid, Molsid logo, SmartID® and ResiScreen® are trademarks of Molsid.

All other trademarks are the property of their respective owners. © 2026 Molsid. All rights reserved.

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