

ResiScreen® Staph



BlaZ β -lactamase activity detection for *Staphylococcus* spp. clinical isolates — Simple, UV-Read

One Assay for *Staphylococcus* spp.

- **Binary, UV readout** — positive or negative
- Covers both ***S. aureus*** and **CoNS***
- Powered by Molsid's SmartID® fluorogenic biotracer

<i>S.a. blaZ+</i>		<i>S.a. blaZ-</i>		
				UV light
				white light
control	test	control	test	



Validated Performance




- **100% overall agreement** with the genotypic reference method **
- Panel representative of clinical strain diversity (n=301)***
- Activity detected across all BlaZ enzyme variants



Fits Your Routine Workflow

- **Minimal hands-on time**
- Reads with a standard 365nm UV light — no extra instrument required
- Runs alongside your routine workflow for AST



* CoNS: coagulase-negative staphylococci

** 95% Clopper-Pearson CI: 98.8–100%

***See Instructions For Use